## Programat S2

### FAQ

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| Program              | Name and Description   | Duration                                 |
|----------------------|--|--|
| <b>1</b><br>Standard | IPS e.max ZirCAD standard sintering of IPS e.max ZirCAD Prime/ MT Multiunit/TMO crowns/bridges up to 14 units and Infrastore, dried restorations | ✓<br>Programat® Dosto Tray<br>9 h 50 min |
| <b>2</b><br>Speed    | IPS e.max ZirCAD Speed – 3 units quick sintering of IPS e.max ZirCAD Prime/ MT Multiunit/TMO crowns/bridges up to 3 units                        | 4 h 25 min                               |
| <b>3</b><br>Speed    | IPS e.max ZirCAD Prime Speed – crown quick sintering of IPS e.max ZirCAD Prime crowns  | 2 h 26 min                               |
| <b>4</b><br>Speed    | IPS e.max ZirCAD MT Multi Speed – 3 units quick sintering of IPS e.max ZirCAD MT Multi crowns/bridges up to 3 units                              | 4 h 25 min                               |
| <b>5</b><br>Speed    | IPS e.max ZirCAD MT Speed – 3 units quick sintering of IPS e.max ZirCAD MT crowns/bridges up to 3 units  | 2 h 30 min                               |
| <b>6</b><br>Speed    | IPS e.max ZirCAD LT Speed – crown quick sintering of IPS e.max ZirCAD LT crowns  | 2 h 55 min                               |
| <b>7</b><br>Speed    | IPS e.max ZirCAD LT Speed – 14 units quick sintering of IPS e.max ZirCAD LT bridges with up to 14 units  | 4 h 30 min                               |
| <b>8</b><br>Speed    | IPS e.max ZirCAD MO Speed – crown quick sintering of IPS e.max ZirCAD MO crown frameworks  | 75 min                                   |
| <b>9</b><br>Speed    | IPS e.max ZirCAD MO Speed – 4 units quick sintering of IPS e.max ZirCAD MO bridge frameworks with up to 4 units                                  | 90 min                                   |
| <b>10</b><br>Speed   | IPS e.max ZirCAD MO Speed – 14 units quick sintering of IPS e.max ZirCAD MO bridge frameworks with up to 14 units                                | 2 h 45 min                               |

|                       |   |                            |            |
|-----------------------|---|----------------------------|------------|
| <b>11</b><br>Standard | IPS e.max ZirCAD MO standard sintering of IPS e.max ZirCAD MO frameworks                                | ✓<br>Programat® Dosto Tray | 7 h 20 min |
| <b>12</b>             | Regeneration firing Program for the regeneration of IPS e.max ZirCAD frameworks after excessive working |                            | 60 min     |
| <b>13</b><br>Standard | Zenostar T standard sintering of Zenostar T   | ✓<br>Programat® Dosto Tray | 9 h 50 min |
| <b>14</b><br>Speed    | Zenostar T crown quick sintering of Zenostar T bridges with up to 14 units                              |                            | 2 h 55 min |
| <b>15</b><br>Speed    | Zenostar T bridge – 14 units quick sintering of Zenostar T bridges with up to 14 units                  |                            | 4 h 30 min |
| <b>16</b><br>Speed    | Zenostar MO crown quick sintering of Zenostar MO crown frameworks                                       |                            | 75 min     |
| <b>17</b><br>Speed    | Zenostar MO bridge – 4 units quick sintering of Zenostar MO bridge frameworks with up to 4 units        |                            | 90 min     |
| <b>18</b><br>Speed    | Zenostar MO bridge – 14 units quick sintering of Zenostar MO bridge frameworks with up to 14 units      |                            | 2 h 45 min |
| <b>19</b><br>Standard | Zenostar MO standard sintering of Zenostar MO frameworks  | ✓<br>Programat® Dosto Tray | 7 h 20 min |

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